

Product datasheet for **RC204168A1V**

Human HOMER3 (NM_004838) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human HOMER3 (NM_004838) AAV Particle
Tag:	Myc-DDK
Symbol:	HOMER3
Synonyms:	HOMER-3; VESL3
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC204168 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGTCCACAGCCAGGGAGCAGCCAATCTTCAGCACACGGGCGCACGTGTTCCAAATTGACCCAGCCACCA
AGCGAACTGGATCCCAGCGGGCAAGCACGCACTCACTGTCTCTATTCTACGATGCCACCCGCAATGT
GTACCGCATCATCAGCATCGGAGGCGCCAAGGCCATCATCAACAGCACTGTCACCTCCCAACATGACCTTC
ACCAAACTTCCAGAAAGTTCGGGCGAGTGGGCCGACAGTCGCGCCAACACAGTCTACGGCTTGGGCTTTG
CCTCTGAACAGCATCTGACACAGTTTGCCGAGAAGTTCAGGAAGTGAAGGAAGCAGCCAGGCTGGCCAG
GGAGAAATCTCAGGATGGCGGGGAGCTCACCAGTCCAGCCCTGGGGCTCGCCTCCCAACAGGTGCCCCCG
AGCCCTCTCGTCAGTGCCAACGGCCCCGCGAGGAAAACTGTTCCGAGCCAGAGCGCTGATGCCCCCG
GCCCCACAGAGCGCGAGCGGCTAAAGAAGATGTTGTCTGAGGGCTCCGTGGGCGAGGTACAGTGGGAGGC
CGAGTTTTTCGCACTGCAGGACAGCAACAACAAGCTGGCAGGCGCCCTGCGAGAGGCCAACGCCCGCGCA
GCCCAGTGGAGGCAGCAGCTGGAGGCTCAGCGTGCAGAGGCCGAGCGGCTGCGGCAGCGGGTGGTGAGC
TGGAGGCTCAGGCAGCTTCAGAGGTGACCCCCACCGGTGAGAAGGAGGGGCTGGGCCAGGGCCAGTCGCT
GGAACAGCTGGAAGCTCTGGTGAAACCAAGGACCAGGAGATTCAGACCCTGAAGAGTCAGACTGGGGG
CCCCGCGAGGCCCTGGAGGCTGCCGAGCGTGAGGAGACTCAGCAGAAGGTGCAGGACCTGGAGACCGCA
ATGCGGAGTTGGAGCACCAGCTGCGGGCGATGGAGCGCAGCCTGGAGGAGGCACGGGCAGAGCGGGAGCG
GGCGCGGGCTGAGGTGGGCCGGCAGCGCAGCTGCTGGACGTCAGGCTGTTTGAGCTGAGTGAGCTGCGT
GAGGGCTGGCCCGCCTGGCTGAGGCTGCGCCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC204168 protein sequence Red=Cloning site Green=Tags(s) MSTAREQPIFSTRAHVFQIDPATKRNWIPAGKHALTVSYFYDATRNVYRIISIGGAKAIINSTVTPNMTF TKTSQKFGQWADSRANTVYGLGFASEQHLTQFAEKFAQEVKEAARLAREKSQDGGELTSPALGLASHQVPP SPLVSANGPGEEKLFRSQSADAPGPTERERLKKMLSEGSVGEVQWEAEFFALQDSNNKLAGALREANAAA AQWRQQLAEQRAEAERLRQVAEAEQAASEVTPTGEKEGLGQGSLEQLEALVQTKDQEIQTLKSQTGG PREALEAAEREETQQKVQDLETRNALEHQLRAMERSLEEARAERERARAQVGRAAQLLDVRLFELSELRL EGLARLAEAAP TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_004838
ORF Size:	1083 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_004838.2 , NP_004829.2
RefSeq Size:	1587 bp
RefSeq ORF:	1086 bp
Locus ID:	9454
UniProt ID:	Q9NSC5
Cytogenetics:	19p13.11
MW:	39.9 kDa