

Product datasheet for **RC204174A1V**

Human SPAG7 (NM_004890) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human SPAG7 (NM_004890) AAV Particle
Tag: Myc-DDK
Symbol: SPAG7
Synonyms: ACRP; FSA-1
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC204174 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGCGGACCTACTGGGCTCCATCCTGAGCTCCATGGAGAAGCCACCCAGCCTCGGTGACCAGGAGACTC
GGCGCAAGGCCCGAGAACAGGCCGCCCGCCTGAAGAACTACAAGAGCAAGAGAAACAACAGAAAGTGGA
GTTTCGTAAAAGGATGGAGAAGGAGGTGTCAGATTTTCATTCAAGACAGTGGGCAGATCAAGAAAAAGTTT
CAGCCAATGAACAAGATCGAGAGGAGCATACTACATGATGTGGTGAAGTGGCTGGCCTGACATCCTTCT
CCTTTGGGAAGATGATGACTGTCGCTATGTGTCATGATCTTCAAAAAGGAGTTTGCACCTCAGATGAAGA
GCTAGACTCTTACCGTCGTGGAGAGGAATGGGACCCCGAGAGGCTGAGGAGAAGCGGAAGCTGAAGGAG
CTGGCCAGAGGCAAGAGGAGGAGGAGGCCAGCCAGCAGGGGCTGTGGTGGTGAAGCCTGCCAGCGACTACA
AGGACAAGTACAGCCACCTCATCGGCAAGGGAGCAGCCAAAGACGCAGCCACATGCTACAGGCCAATAA
GACCTACGGCTGTGTGCCGTGGCCAATAAGAGGGACACACGCTCCATTGAAGAGGCTATGAATGAGATC
AGAGCCAAGAAGCGTCTGCGGCAGAGTGGGAAGAGTTGCCGCCAACCTCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC204174 protein sequence
Red=Cloning site Green=Tags(s)

MADLLGSILSSMEKPPSLGDQETRRKAREQAARLKKLQEKEQKQKVEFRKRMEKEVSDFIQDSGQIKKKF
QPMNKIERSILHDVVEVAGLTSFSFGEDDDCRYVMIFFKKEFAPSDEELDSYRRGEEWDPQKAEEKRKLKE
LAQRQEEEAQQGPVVVSPADYKDKYSHLIGKAAKDAHMLQANKTYGCVPVANKRDRSIEEAMNEI
RAKKRLRQSGEELPPTS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV



[View online »](#)

Species:	Human
Serotype:	AAV-2
ACCN:	NM_004890
ORF Size:	681 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_004890.2
RefSeq Size:	1038 bp
RefSeq ORF:	684 bp
Locus ID:	9552
UniProt ID:	O75391
Cytogenetics:	17p13.2
MW:	26 kDa