

Product datasheet for **RC204178L2V**

GNB2 (NM_005273) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	GNB2 (NM_005273) Human Tagged ORF Clone Lentiviral Particle
Symbol:	GNB2
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_005273
ORF Size:	1020 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204178).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_005273.2
RefSeq Size:	1681 bp
RefSeq ORF:	1023 bp
Locus ID:	2783
UniProt ID:	P62879
Cytogenetics:	7q22.1
Domains:	WD40
Protein Pathways:	Chemokine signaling pathway
MW:	37.3 kDa



[View online »](#)

Gene Summary:

Heterotrimeric guanine nucleotide-binding proteins (G proteins), which integrate signals between receptors and effector proteins, are composed of an alpha, a beta, and a gamma subunit. These subunits are encoded by families of related genes. This gene encodes a beta subunit. Beta subunits are important regulators of alpha subunits, as well as of certain signal transduction receptors and effectors. This gene contains a trinucleotide (CCG) repeat length polymorphism in its 5' UTR. [provided by RefSeq, Jul 2008]