

Product datasheet for **RC204184L1V**

NFIC (NM_005597) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	NFIC (NM_005597) Human Tagged ORF Clone Lentiviral Particle
Symbol:	NFIC
Synonyms:	CTF; CTF5; NF-I; NFI
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_005597
ORF Size:	1284 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204184).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_005597.2
RefSeq Size:	7823 bp
RefSeq ORF:	1287 bp
Locus ID:	4782
UniProt ID:	P08651
Cytogenetics:	19p13.3
Protein Families:	Transcription Factors
MW:	47.6 kDa



[View online »](#)

Gene Summary:

The protein encoded by this gene belongs to the CTF/NF-I family. These are dimeric DNA-binding proteins, and function as cellular transcription factors and as replication factors for adenovirus DNA replication. Alternatively spliced transcript variants encoding different isoforms have been described for this gene. [provided by RefSeq, Oct 2011]