

Product datasheet for **RC204184L3V**

NFIC (NM_005597) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | NFIC (NM_005597) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | NFIC |
| Synonyms: | CTF; CTF5; NF-I; NFI |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-Myc-DDK-P2A-Puro (PS100092) |
| Tag: | Myc-DDK |
| ACCN: | NM_005597 |
| ORF Size: | 1284 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC204184). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_005597.2 |
| RefSeq Size: | 7823 bp |
| RefSeq ORF: | 1287 bp |
| Locus ID: | 4782 |
| UniProt ID: | P08651 |
| Cytogenetics: | 19p13.3 |
| Protein Families: | Transcription Factors |
| MW: | 47.6 kDa |



[View online »](#)

Gene Summary:

The protein encoded by this gene belongs to the CTF/NF-I family. These are dimeric DNA-binding proteins, and function as cellular transcription factors and as replication factors for adenovirus DNA replication. Alternatively spliced transcript variants encoding different isoforms have been described for this gene. [provided by RefSeq, Oct 2011]