

Product datasheet for **RC204186L2V**

PEDF (SERPINF1) (NM_002615) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | PEDF (SERPINF1) (NM_002615) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | PEDF |
| Synonyms: | EPC-1; OI6; OI12; PEDF; PIG35 |
| Mammalian Cell Selection: | None |
| Vector: | pLenti-C-mGFP (PS100071) |
| Tag: | mGFP |
| ACCN: | NM_002615 |
| ORF Size: | 1254 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC204186). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_002615.4 |
| RefSeq Size: | 1552 bp |
| RefSeq ORF: | 1257 bp |
| Locus ID: | 5176 |
| UniProt ID: | P36955 |
| Cytogenetics: | 17p13.3 |
| Domains: | SERPIN |
| Protein Families: | Druggable Genome, Secreted Protein |



[View online »](#)

MW: 46.4 kDa

Gene Summary: This gene encodes a member of the serpin family that does not display the serine protease inhibitory activity shown by many of the other serpin proteins. The encoded protein is secreted and strongly inhibits angiogenesis. In addition, this protein is a neurotrophic factor involved in neuronal differentiation in retinoblastoma cells. Mutations in this gene were found in individuals with osteogenesis imperfecta, type VI. [provided by RefSeq, Aug 2016]