

## Product datasheet for **RC204204L2V**

### SIL1 (NM\_001037633) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	SIL1 (NM_001037633) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SIL1
Synonyms:	BAP; MSS; ULG5
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_001037633
ORF Size:	1383 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204204).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001037633.1</a> , <a href="#">NP_001032722.1</a>
RefSeq Size:	1996 bp
RefSeq ORF:	1386 bp
Locus ID:	64374
UniProt ID:	<a href="#">Q9H173</a>
Cytogenetics:	5q31.2
Protein Families:	Protease, Secreted Protein, Transmembrane
MW:	52.1 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a resident endoplasmic reticulum (ER), N-linked glycoprotein with an N-terminal ER targeting sequence, 2 putative N-glycosylation sites, and a C-terminal ER retention signal. This protein functions as a nucleotide exchange factor for another unfolded protein response protein. Mutations in this gene have been associated with Marinesco-Sjogren syndrome. Alternate transcriptional splice variants have been characterized. [provided by RefSeq, Jul 2008]