

## Product datasheet for **RC204208L3V**

### Dysbindin (DTNBP1) (NM\_032122) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Dysbindin (DTNBP1) (NM_032122) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Dysbindin
Synonyms:	BLOC1S8; DBND; HPS7; My031; SDY
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_032122
ORF Size:	1053 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204208).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_032122.3</a>
RefSeq Size:	1429 bp
RefSeq ORF:	1056 bp
Locus ID:	84062
UniProt ID:	<a href="#">Q96EV8</a>
Cytogenetics:	6p22.3
Domains:	DUF546
Protein Families:	Druggable Genome


[View online »](#)

**MW:** 39.5 kDa

**Gene Summary:** This gene encodes a protein that may play a role in organelle biogenesis associated with melanosomes, platelet dense granules, and lysosomes. A similar protein in mouse is a component of a protein complex termed biogenesis of lysosome-related organelles complex 1 (BLOC-1), and binds to alpha- and beta-dystrobrevins, which are components of the dystrophin-associated protein complex (DPC). Mutations in this gene are associated with Hermansky-Pudlak syndrome type 7. This gene may also be associated with schizophrenia. Multiple transcript variants encoding distinct isoforms have been identified for this gene. [provided by RefSeq, Jul 2008]