

Product datasheet for **RC204210L4V**

TNIP1 (NM_006058) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	TNIP1 (NM_006058) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TNIP1
Synonyms:	ABIN-1; NAF1; nip40-1; VAN
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_006058
ORF Size:	1908 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204210).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_006058.3
RefSeq Size:	2933 bp
RefSeq ORF:	1911 bp
Locus ID:	10318
UniProt ID:	Q15025
Cytogenetics:	5q33.1
MW:	71.9 kDa



[View online »](#)

Gene Summary:

This gene encodes an A20-binding protein which plays a role in autoimmunity and tissue homeostasis through the regulation of nuclear factor kappa-B activation. Mutations in this gene have been associated with psoriatic arthritis, rheumatoid arthritis, and systemic lupus erythematosus. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Nov 2011]