

Product datasheet for **RC204220L3V**

NOXO1 (NM_144603) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	NOXO1 (NM_144603) Human Tagged ORF Clone Lentiviral Particle
Symbol:	NOXO1
Synonyms:	P41NOX; P41NOXA; P41NOXB; P41NOXC; SH3PXD5; SNX28
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_144603
ORF Size:	1110 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204220).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_144603.2
RefSeq Size:	1636 bp
RefSeq ORF:	1113 bp
Locus ID:	124056
UniProt ID:	Q8NFA2
Cytogenetics:	16p13.3
MW:	40.7 kDa



[View online »](#)

Gene Summary:

This gene encodes an NADPH oxidase (NOX) organizer, which positively regulates NOX1 and NOX3. The protein contains a PX domain and two SH3 domains. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Jun 2012]