

Product datasheet for **RC204238A1V**

Human Glutamine Synthetase (GLUL) (NM_001033044) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human Glutamine Synthetase (GLUL) (NM_001033044) AAV Particle
Tag:	Myc-DDK
Symbol:	Glutamine Synthetase
Synonyms:	GLNS; GS; PIG43; PIG59
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC204238 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGACCACCTCAGCAAGTTCCCACTTAAATAAAGGCATCAAGCAGGTGTACATGTCCCTGCCTCAGGGTG
AGAAAGTCCAGGCCATGTATATCTGGATCGATGGTACTGGAGAAGGACTGCGCTGCAAGACCCGGACCCT
GGACAGTGAGCCCAAGTGTGTGAAGAGTTGCCCTGAGTGGAAATTCGATGGCTCTAGTACTTTACAGTCT
GAGGGTTCCAACAGTGACATGTATCTCGTGCCTGCTGCCATGTTTCGGGACCCCTCCGTAAGGACCCTA
ACAAGCTGGTGTATGTGAAGTTTCAAGTACAATCGAAGGCCTGCAGAGACCAATTTGAGGCACACCTG
TAAACGGATAATGGACATGGTGAGCAACCAGCACCCCTGGTTTGGCATGGAGCAGGAGTATACCCCTCATG
GGGACAGATGGGCACCCCTTTGGTTGGCCTTCCAACGGCTTCCAGGGCCCCAGGGTCCATATTACTGTG
GTGTGGGAGCAGACAGAGCCTATGGCAGGGACATCGTGGAGGCCATTACCGGGCCTGCTTGTATGCTGG
AGTCAAGATTGCGGGGACTAATGCCGAGGTGATGCCTGCCAGTGGGAATTTAGATTGGACCTTGTGAA
GGAATCAGCATGGGAGATCATCTCTGGTGGCCGTTTCATCTGCATCGTGTGTGAAGACTTTGGAG
TGATAGCAACCTTTGATCCTAAGCCATTCTGGGAAGTGAATGGTGCAGGCTGCCATACCAACTTCAG
CACCAAGGCCATGCGGGAGGAGAATGGTCTGAAGTACATCGAGGAGGCCATTGAGAACTAAGCAAGCGG
CACCAGTACCACATCCGTGCCTATGATCCCAAGGAGGCCTGGACAATGCCGACGTCTAACTGGATTCC
ATGAAACCTCCAACATCAACGACTTTTCTGCTGGTGTAGCCAATCGTAGCGCCAGCATACGATTCCCCG
GACTGTTGGCCAGGAGAAGAAGGTTACTTTGAAGATCGTCGCCCTCTGCCAACTGCGACCCCTTTTCG
GTGACAGAAGCCCTCATCCGCACGTGTCTTCTCAATGAAACCGGCGATGAGCCCTTCAGTACAAAAAT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC204238 protein sequence Red=Cloning site Green=Tags(s) MTTSASSHLNKGIKQVYMSLPQGEKVQAMYIWIDGTGEGLRCKTRTLTLDSEPKCVEELPEWNFDGSSTLQS EGNSNDMYLVPAAMFRDPFRKDPNKLVLCEVFKYNRRPAETNLRHTCKRIMDMVSNQHPWFGMEQEYTLT GTDGHPFGWPSNGFPGPQGPYYCGVGADRAYGRDIVEAHYRACLYAGVKIAGTNAEVMPAQWFEFIGPCE GISMGDHLWVARFILHRVCEDFGVIAITFDPKPIPGNWNAGAGCHTNFSTKAMREENGLKYIEEAIEKLSKR HQYHIRAYDPKGGLDNARRLTGFHETSINIDFSAGVANRSASIRIPRTVGQEKKG YFEDRRPSANCDPFS VTEALIRTCLLNETGDEPFQYKN TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_001033044
ORF Size:	1119 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_001033044.1
RefSeq Size:	4381 bp
RefSeq ORF:	1122 bp
Locus ID:	2752
UniProt ID:	P15104
Cytogenetics:	1q25.3
MW:	42.1 kDa

Product images: