

Product datasheet for **RC204241L4V**

RRP36 (NM_033112) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	RRP36 (NM_033112) Human Tagged ORF Clone Lentiviral Particle
Symbol:	RRP36
Synonyms:	C6orf153; dj20C7.4
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_033112
ORF Size:	777 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204241).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_033112.2
RefSeq Size:	1177 bp
RefSeq ORF:	780 bp
Locus ID:	88745
UniProt ID:	Q96EU6
Cytogenetics:	6p21.1
MW:	30.3 kDa



[View online »](#)

Gene Summary:

RRP36 functions at an early stage in the processing of 35S preribosomal RNA into the mature 18S species (Gerus et al., 2010 [PubMed 20038530]).[supplied by OMIM, Jul 2010]