

Product datasheet for **RC204243A1V**

Human Tcl1 (TCL1A) (NM_021966) AAV Particle

Product data:

Product Type: AAV Particles

Product Name: Human Tcl1 (TCL1A) (NM_021966) AAV Particle

Tag: Myc-DDK

Symbol: Tcl1

Synonyms: TCL1

Mammalian Cell Selection: None

Vector: pAAV-AC-Myc-DDK (PS100089)

ORF Nucleotide Sequence: >RC204243 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGCCGAGTGCCCGACACTCGGGGAGGCAGTCACCGACACCCGGACCGCCTGTGGGCCTGGGAGAAGT
TCGTGTATTTGGACGAGAAGCAGCAGCCTGGCTGCCCTTAACCATCGAGATAAAGGATAGGTTACAGTT
ACGGGTGCTCTTGGCTCGGGAAGACGTCGTCTGGGGAGGCCTATGACCCCAACCCAGATAGGCCCAAGC
CTGCTGCCTATCATGTGGCAGCTCTACCCTGATGGACGATACCGATCCTCAGACTCCAGTTTCTGGCGCT
TAGTGTACCACATCAAGATTGACGGCGTGGAGGACATGCTTCTCGAGCTGCTGCCAGATGAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC204243 protein sequence
Red=Cloning site Green=Tags(s)

MAECPTLGEAVTDHPDRLWAWKFFVYLDEKQHAWLPLTIEIKDRLQLRVLLRREDVVLGRPMTPTQIGPS
LLPIMWQLYPDGRYSSDSSFWRVYHIKIDGVEDMLLELLPDD

TRTRPLE**QKL**ISEEDLAANDILDYKDDDDK**V**

Species: Human

Serotype: AAV-2

ACCN: NM_021966

ORF Size: 342 bp



[View online »](#)

Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_021966.1
RefSeq Size:	1410 bp
RefSeq ORF:	345 bp
Locus ID:	8115
UniProt ID:	P56279
Cytogenetics:	14q32.13
MW:	13.5 kDa