

Product datasheet for **RC204244L3V**

SEPT6 (SEPTIN6) (NM_145802) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	SEPT6 (SEPTIN6) (NM_145802) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SEPTIN6
Synonyms:	SEP2; SEPT2; SEPT6
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_145802
ORF Size:	1287 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204244).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_145802.2
RefSeq Size:	2564 bp
RefSeq ORF:	1290 bp
Locus ID:	23157
UniProt ID:	Q14141
Cytogenetics:	Xq24
Domains:	GTP_CDC
MW:	49.2 kDa



[View online »](#)

Gene Summary:

This gene is a member of the septin family of GTPases. Members of this family are required for cytokinesis. One version of pediatric acute myeloid leukemia is the result of a reciprocal translocation between chromosomes 11 and X, with the breakpoint associated with the genes encoding the mixed-lineage leukemia and septin 2 proteins. This gene encodes four transcript variants encoding three distinct isoforms. An additional transcript variant has been identified, but its biological validity has not been determined. [provided by RefSeq, Jul 2008]