

Product datasheet for **RC204249L4V**

NFX1 (NM_002504) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	NFX1 (NM_002504) Human Tagged ORF Clone Lentiviral Particle
Symbol:	NFX1
Synonyms:	NFX2; TEG-42; Tex42
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_002504
ORF Size:	3360 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204249).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_002504.3
RefSeq Size:	4721 bp
RefSeq ORF:	3363 bp
Locus ID:	4799
UniProt ID:	Q12986
Cytogenetics:	9p13.3
Domains:	zf-NF-X1, R3H
Protein Families:	Druggable Genome, Transcription Factors



[View online »](#)

MW: 124.4 kDa

Gene Summary: MHC class II gene expression is controlled primarily at the transcriptional level by transcription factors that bind to the X and Y boxes, two highly conserved elements in the proximal promoter of MHC class II genes. The protein encoded by this gene is a transcriptional repressor capable of binding to the conserved X box motif of HLA-DRA and other MHC class II genes in vitro. The protein may play a role in regulating the duration of an inflammatory response by limiting the period in which class II MHC molecules are induced by IFN-gamma. Three alternative splice variants, each of which encodes a different isoform, have been identified. [provided by RefSeq, Jul 2008]