

Product datasheet for **RC204251**

RPL28 (NM_000991) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: RPL28 (NM_000991) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: RPL28
Synonyms: L28
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >RC204251 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTCTGCGCATCTGCAATGGATGGTCTGTCGGAAGTCTCCAGTTTCCTGATCAAGAGGAATAAGCAGA
CCTACAGCACTGAGCCCAATAACTTGAAGGCCCGCAATTCCTTCCGCTACAACGGACTGATTCACCGCAA
GACTGTGGCGTGGAGCCGGCAGCCGACGGCAAAGGTGTCGTGGTGGTCATTAAGCGGAGATCCGGCCAG
CGGAAGCCTGCCACCTCCTATGTGCGGACCACCATCAACAAGAATGCTCGCGCCACGCTCAGCAGCATCA
GACACATGATCCGCAAGAACAAGTACCGCCCCGACCTGCGCATGGCAGCCATCCGCAGGGCCAGCGCCAT
CCTGCGCAGCCAGAAGCCTGTGATGGTGAAGAGGAAGCGGACCCGCCCAAGAGCTCC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC204251 protein sequence
Red=Cloning site Green=Tags(s)

MSAHLQWMVVRNCSSFLIKRNKQTYSTEPNNLKARNSFRYNGLIHRKTVGVPEAADGKGVVVVVIKRRSQ
RKPATSYVRTTINKNARATLSSIRHMIRKNKYRPDLRMAAIRRASAILRSQKPVVMVKRKRTRPTKSS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk6059_e12.zip

Restriction Sites: SgfI-MluI



[View online »](#)

Cloning Scheme:


ACCN: NM_000991

ORF Size: 411 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_000991.5](#)

RefSeq Size: 4245 bp

RefSeq ORF: 414 bp

Locus ID: 6158

UniProt ID: [P46779](#)

Cytogenetics: 19q13.42


Domains: Ribosomal_L28e

Protein Pathways: Ribosome


MW: 15.7 kDa

Gene Summary: Ribosomes, the organelles that catalyze protein synthesis, consist of a small 40S subunit and a large 60S subunit. Together these subunits are composed of 4 RNA species and approximately 80 structurally distinct proteins. This gene encodes a ribosomal protein that is a component of the 60S subunit. The protein belongs to the L28E family of ribosomal proteins. It is located in the cytoplasm. Variable expression of this gene in colorectal cancers compared to adjacent normal tissues has been observed, although no correlation between the level of expression and the severity of the disease has been found. As is typical for genes encoding ribosomal proteins, there are multiple processed pseudogenes of this gene dispersed through the genome. Alternative splicing results in multiple transcript variants encoding distinct isoforms.[provided by RefSeq, Oct 2008]

Product images:



Circular map for RC204251



Western blot validation of overexpression lysate (Cat# [LY400358]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC204251 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).