

Product datasheet for **RC204266L2V**

FAM120B (NM_032448) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	FAM120B (NM_032448) Human Tagged ORF Clone Lentiviral Particle
Symbol:	FAM120B
Synonyms:	CCPG; dj894D12.1; KIAA1838; PGCC1; SAN1
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_032448
ORF Size:	2730 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204266).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_032448.1 , NP_115824.1
RefSeq Size:	5211 bp
RefSeq ORF:	2733 bp
Locus ID:	84498
UniProt ID:	Q96EK7
Cytogenetics:	6q27
MW:	103.8 kDa



[View online »](#)

Gene Summary:

Functions as a transactivator of PPARG and ESR1. Functions in adipogenesis through PPARG activation (By similarity).[UniProtKB/Swiss-Prot Function]