

Product datasheet for **RC204276L4V**

ZNF307 (ZKSCAN4) (NM_019110) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ZNF307 (ZKSCAN4) (NM_019110) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ZNF307
Synonyms:	P1P373C6; ZNF307; ZNF427; ZSCAN36
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_019110
ORF Size:	1635 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204276).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_019110.3
RefSeq Size:	2307 bp
RefSeq ORF:	1638 bp
Locus ID:	387032
UniProt ID:	Q969J2
Cytogenetics:	6p22.1
Domains:	KRAB, LER, zf-C2H2
Protein Families:	Transcription Factors



[View online »](#)

MW: 61.6 kDa

Gene Summary: May be involved in the transcriptional activation of MDM2 and EP300 genes.
[UniProtKB/Swiss-Prot Function]