

## Product datasheet for **RC204300L2V**

### WDR21A (DCAF4) (NM\_015604) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	WDR21A (DCAF4) (NM_015604) Human Tagged ORF Clone Lentiviral Particle
Symbol:	WDR21A
Synonyms:	WDR21; WDR21A
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_015604
ORF Size:	1485 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204300).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_015604.2</a>
RefSeq Size:	2563 bp
RefSeq ORF:	1488 bp
Locus ID:	26094
UniProt ID:	<a href="#">Q8WV16</a>
Cytogenetics:	14q24.2
Domains:	WD40
MW:	55.5 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a WD repeat-containing protein that interacts with the Cul4-Ddb1 E3 ligase macromolecular complex. Multiple alternatively spliced transcript variants have been found for this gene. [provided by RefSeq, Jul 2009]