

Product datasheet for **RC204313L3V**

PRR13 (NM_018457) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	PRR13 (NM_018457) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PRR13
Synonyms:	TXR1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_018457
ORF Size:	444 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204313).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_018457.2
RefSeq Size:	1262 bp
RefSeq ORF:	447 bp
Locus ID:	54458
UniProt ID:	Q9NZ81
Cytogenetics:	12q13.13
MW:	15.4 kDa



[View online »](#)

Gene Summary:

Negatively regulates TSP1 expression at the level of transcription. This down-regulation was shown to reduce taxane-induced apoptosis.[UniProtKB/Swiss-Prot Function]