

## Product datasheet for **RC204316L4V**

### SLC35C2 (NM\_173073) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	SLC35C2 (NM_173073) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SLC35C2
Synonyms:	BA394O2.1; C20orf5; CGI-15; OVCOV1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_173073
ORF Size:	1032 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204316).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_173073.2</a>
RefSeq Size:	2262 bp
RefSeq ORF:	1035 bp
Locus ID:	51006
UniProt ID:	<a href="#">Q9NQQ7</a>
Cytogenetics:	20q13.12
Protein Families:	Druggable Genome, Transmembrane
MW:	38.2 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a member of the triose-phosphate transporter protein family. This gene is regulated by oxygen tension, is induced in hypoxic trophoblast cells, and is overexpressed in ovarian cancer. Alternative splicing results in multiple transcript variants. A pseudogene of this gene has been defined on the X chromosome. [provided by RefSeq, Jul 2013]