

## Product datasheet for **RC204317L1V**

### GLUT8 (SLC2A8) (NM\_014580) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	GLUT8 (SLC2A8) (NM_014580) Human Tagged ORF Clone Lentiviral Particle
Symbol:	GLUT8
Synonyms:	GLUT8; GLUTX1
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_014580
ORF Size:	1431 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204317).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_014580.3</a>
RefSeq Size:	2176 bp
RefSeq ORF:	1434 bp
Locus ID:	29988
UniProt ID:	<a href="#">Q9NY64</a>
Cytogenetics:	9q33.3
Protein Families:	Transmembrane
MW:	50.8 kDa



[View online »](#)

**Gene Summary:**

This gene belongs to the solute carrier 2A family, which includes intracellular glucose transporters. Based on sequence comparison, the glucose transporters are grouped into three classes and this gene is a member of class II. The encoded protein, like other members of the family, contains several conserved residues and motifs and 12 transmembrane domains with both amino and carboxyl ends being on the cytosolic side of the membrane. Alternatively spliced transcript variants have been described for this gene. [provided by RefSeq, Nov 2012]