

### Product datasheet for RC204325L1

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## beta glucuronidase (GUSB) (NM\_000181) Human Tagged Lenti ORF Clone

**Product data:** 

**Product Type:** Expression Plasmids

**Product Name:** beta glucuronidase (GUSB) (NM\_000181) Human Tagged Lenti ORF Clone

Tag: Myc-DDK

Symbol: beta glucuronidase

**Synonyms:** BG; MPS7

Mammalian Cell None

Selection:

Vector:pLenti-C-Myc-DDK (PS100064)E. coli Selection:Chloramphenicol (34 ug/mL)


ORF Nucleotide The ORF insert of this clone is exactly the same as(RC204325).

Sequence:

**Restriction Sites:** Sgfl-Mlul

**Cloning Scheme:** 





<sup>\*</sup> The last codon before the Stop codon of the ORF.

ACCN: NM\_000181

ORF Size: 1953 bp



#### beta glucuronidase (GUSB) (NM\_000181) Human Tagged Lenti ORF Clone - RC204325L1

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:** 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

**RefSeq:** <u>NM 000181.2</u>

 RefSeq Size:
 2321 bp

 RefSeq ORF:
 1956 bp

 Locus ID:
 2990

 UniProt ID:
 P08236

Cytogenetics: 7q11.21

**Domains:** Glyco\_hydro\_2, Glyco\_hydro\_2\_C, Glyco\_hydro\_2\_N

**Protein Pathways:** Drug metabolism - other enzymes, Glycosaminoglycan degradation, Lysosome, Metabolic

pathways, Pentose and glucuronate interconversions, Porphyrin and chlorophyll metabolism,

Starch and sucrose metabolism

Druggable Genome, Transmembrane

**MW:** 74.7 kDa

**Protein Families:** 


**Gene Summary:** This gene encodes a hydrolase that degrades glycosaminoglycans, including heparan sulfate,

dermatan sulfate, and chondroitin-4,6-sulfate. The enzyme forms a homotetramer that is localized to the lysosome. Mutations in this gene result in mucopolysaccharidosis type VII. Alternative splicing results in multiple transcript variants. There are many pseudogenes of this


locus in the human genome.[provided by RefSeq, May 2014]



# **Product images:**



Circular map for RC204325L1



Double digestion of RC204325L1 using Sgfl and Mlul