

## Product datasheet for **RC204332L2V**

### **SMPDL3A (NM\_006714) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	SMPDL3A (NM_006714) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SMPDL3A
Synonyms:	ASM3A; ASML3a; yR36GH4.1
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_006714
ORF Size:	1359 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204332).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_006714.2</a>
RefSeq Size:	1987 bp
RefSeq ORF:	1362 bp
Locus ID:	10924
UniProt ID:	<a href="#">Q92484</a>
Cytogenetics:	6q22.31
Domains:	Metallophos
Protein Families:	Secreted Protein



[View online »](#)

**MW:** 51.3 kDa

**Gene Summary:** Has in vitro nucleotide phosphodiesterase activity with nucleoside triphosphates, such as ATP (PubMed:25288789, PubMed:26783088). Has in vitro activity with p-nitrophenyl-TMP (PubMed:25288789). Has lower activity with nucleoside diphosphates, and no activity with nucleoside monophosphates (PubMed:25288789, PubMed:26783088). Has in vitro activity with CDP-choline, giving rise to CMP and phosphocholine. Has in vitro activity with CDP-ethanolamine (PubMed:26783088). Does not have sphingomyelin phosphodiesterase activity (PubMed:25288789, PubMed:26783088).[UniProtKB/Swiss-Prot Function]