

Product datasheet for **RC204338L3V**

AP4B1 (NM_006594) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	AP4B1 (NM_006594) Human Tagged ORF Clone Lentiviral Particle
Symbol:	AP4B1
Synonyms:	BETA-4; CPSQ5; SPG47
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_006594
ORF Size:	2217 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204338).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_006594.1
RefSeq Size:	2791 bp
RefSeq ORF:	2220 bp
Locus ID:	10717
UniProt ID:	Q9Y6B7
Cytogenetics:	1p13.2
Domains:	Adaptin_N
Protein Pathways:	Lysosome



[View online »](#)

MW: 83.3 kDa

Gene Summary: This gene encodes a subunit of a heterotetrameric adapter-like complex 4 that is involved in targeting proteins from the trans-Golgi network to the endosomal-lysosomal system. Mutations in this gene are associated with cerebral palsy spastic quadriplegic type 5 (CPSQ5) disorder. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Dec 2011]