

## Product datasheet for **RC204342L2V**

### ARH3 (ADPRHL2) (NM\_017825) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	ARH3 (ADPRHL2) (NM_017825) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ARH3
Synonyms:	ADPRHL2; ARH3; CONDSIAS
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_017825
ORF Size:	1089 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204342).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_017825.1</a>
RefSeq Size:	1701 bp
RefSeq ORF:	1092 bp
Locus ID:	54936
UniProt ID:	<a href="#">Q9NX46</a>
Cytogenetics:	1p34.3
Domains:	ADP-ribosyl_GH
MW:	38.9 kDa


[View online »](#)

**Gene Summary:**

This gene encodes a member of the ADP-ribosylglycohydrolase family. The encoded enzyme catalyzes the removal of ADP-ribose from ADP-ribosylated proteins. This enzyme localizes to the mitochondria, in addition to the nucleus and cytoplasm.[provided by RefSeq, Feb 2009]