

Product datasheet for **RC204347A1V**

Human Syntaxin 16 (STX16) (NM_003763) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human Syntaxin 16 (STX16) (NM_003763) AAV Particle
Tag:	Myc-DDK
Symbol:	Syntaxin 16
Synonyms:	SYN16
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC204347 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCCACCAGGCGTTTAACCGACGCTTTCTTGTGTGCGGAATAATTCCATCCAAAACCGGCAGCTGT
TAGCCGAGCAACTTGCTGATGACCGTATGGCACTGGTGTGAGGCATCAGCTTAGATCCAGAAGCAGCGAT
TGGTGTGACAAAACGGCCACCTCCTAAGTGGGTGGATGGAGTGGATGAAATTCAGTATGATGTTGGCCGG
ATTAAGCAGAAGATGAAAGAATTGGCCAGCCTTCATGACAAGCATTTAAACAGACCCACCCTGGATGACA
GCAGCGAAGAGGAACATGCCATTGAGATAACTACCCAAGAGATCACTCAGCTCTTCCACAGGTGCCAGCG
TGCCGTGCAGGCCCTGCCGAGCCGGCCCGGGCCTGCTCCGAGCAGGAGGGGCGGCTGCTTGGGAACGTG
GTGGCCTCGCTGGCGCAGGCCCTGCAGGAACCTCCACCAGCTTCCGGCACGCACAGTCAGGCTACCTCA
AACGCATGAAGAATCGAGAGGAAAGATCCCAGCATTTTTTCGACACATCAGTACCACTAATGGATGATGG
AGACGATAACACTCTTTACCATCGGGGTTTACAGAGGACCAGTTAGTTCTGGTGGAGCAGAACACACTG
ATGGTGAAGAGCGGGAACGAGAGATTCGCCAGATTGTACAGTCCATTTCTGACCTGAATGAAATATTCA
GGGACTTAGGGGCGATGATTGTAGAACAGGGTACAGTCCTTGACAGAATTGACTATAACGTTGAACAGTC
CTGTATCAAACTGAAGATGGTTTGAACAGCTTCACAAGGCAGAACAGTATCAAAGAAGAATCGGAAG
ATGCTTGTGATTTAATATTATTGTCATCATCATTGTGCTCATTGTTGCTCCTGCTTGGCGTGAAGTCTC
GA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC204347 protein sequence Red=Cloning site Green=Tags(s) MATRRLTDAFLLLRNNSIQNRQLLAEQLADDRMALVSGISLDPEAAIGVTKRPPPKWVDGVDEIQYDVGR IKQKMKELASLHDKHLNRPTLDDSSSEEEHAIEITTQEITQLFHRCQRAVQALPSRARACSEGEGRLLGNV VASLAQALQELSTSFRAHQSGYLKRMKNREERSQHFFDTSVPLMDDGDDNTLYHRGFTEDQLVLVEQNTL MVEEREREIRQIVQSISDLNEIFRDLGAMIVEQGTVLDRIDYNVEQSCIKTEDGLKQLHKAQYQKKNRK MLVILILFVIIIVLIVLVGVKSR TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_003763
ORF Size:	912 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_003763.3
RefSeq Size:	4915 bp
RefSeq ORF:	915 bp
Locus ID:	8675
UniProt ID:	O14662
Cytogenetics:	20q13.32
MW:	34.8 kDa