

Product datasheet for **RC204348L4V**

CTHRC1 (NM_138455) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	CTHRC1 (NM_138455) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CTHRC1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_138455
ORF Size:	729 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204348).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_138455.2 , NP_612464.1
RefSeq Size:	1279 bp
RefSeq ORF:	732 bp
Locus ID:	115908
UniProt ID:	Q96CG8
Cytogenetics:	8q22.3
Protein Families:	Secreted Protein
MW:	26.2 kDa



[View online »](#)

Gene Summary:

This locus encodes a protein that may play a role in the cellular response to arterial injury through involvement in vascular remodeling. Mutations at this locus have been associated with Barrett esophagus and esophageal adenocarcinoma. Alternatively spliced transcript variants have been described. [provided by RefSeq, Jan 2012]