

Product datasheet for **RC204358L2V**

TCP1 epsilon (CCT5) (NM_012073) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	TCP1 epsilon (CCT5) (NM_012073) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TCP1 epsilon
Synonyms:	CCT-epsilon; CCTE; HEL-S-69; PNAS-102; TCP-1-epsilon
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_012073
ORF Size:	1623 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204358).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_012073.3
RefSeq Size:	3403 bp
RefSeq ORF:	1626 bp
Locus ID:	22948
UniProt ID:	P48643
Cytogenetics:	5p15.2
Domains:	cpn60_TCP1
Protein Families:	Druggable Genome



[View online »](#)

MW: 59.7 kDa

Gene Summary: The protein encoded by this gene is a molecular chaperone that is a member of the chaperonin containing TCP1 complex (CCT), also known as the TCP1 ring complex (TRiC). This complex consists of two identical stacked rings, each containing eight different proteins. Unfolded polypeptides enter the central cavity of the complex and are folded in an ATP-dependent manner. The complex folds various proteins, including actin and tubulin. Mutations in this gene cause hereditary sensory and autonomic neuropathy with spastic paraplegia (HSNSP). Alternative splicing results in multiple transcript variants. Related pseudogenes have been identified on chromosomes 5 and 13. [provided by RefSeq, Apr 2015]