

Product datasheet for **RC204360A1V**

Human HNRPAB (HNRNPAB) (NM_031266) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human HNRPAB (HNRNPAB) (NM_031266) AAV Particle
Tag:	Myc-DDK
Symbol:	HNRPAB
Synonyms:	ABBP1; HNRPAB
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC204360 representing NM_031266 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGC**

ATGTCGGAAGCGGGCGAGGAGCAGCCCATGGAGACGACGGGCGCCACCGAGAACGGACATGAGGCCGTCC
CCGAAGGCGAGTCGCGGCCGGGGCTGGCAGGGCGCCGCGCGGGGGCTGGAGGCGCGACCGCGGCC
CCCGAGCGGGAATCAGAACGGCGCCGAGGGCGACCATCAACGCCAGCAAGAACGAGGAGGACGCGGGA
AAAATGTTCTGTTGGTGGCCTGAGCTGGGATACTAGCAAAAAAGATTTAAAGACTATTTTACTAAATTTG
GAGAGGTCGTTGACTGTACAATAAAATGGATCCCAACTGGACGGTCAAGAGGTTTGGGTTTATCCT
GTTCAAAGATGCAGCCAGTGTGGAGAAGGTCCTAGACCAGAAGGAGCACAGGCTGGATGGCCGTGTCATT
GACCCTAAAAAGGCCATGGCTATGAAGAAGGACCCGGTGAAGAAAATCTTCGTTGGGGGTCTGAATCCTG
AAGCCACTGAGGAAAAGATCAGGGAGTACTTTGGCGAGTTTGGGAGATTGAGGCCATTGAATTGCCAAT
GGATCCAAAGTTGAACAAAAGACGAGGTTTGTGTTTATCACCTTTAAAGAAGAAGAACCCGTGAAGAAG
GTTCTGGAGAAAAAGTTCCATACTGTCACTGGAAGCAAGTGTGAGATCAAGGTGGCCAGCCAAAGAAG
TCTATCAGCAGCAGCAGTATGGCTCTGGGGCCGTGGAACCGCAACCGAGGGAACCGAGGCAGCGGAGG
TGGTGGTGGAGGTGAGGTGAGTCAAGTTCAGGCTACGGGCTACGGCAACTACTGGAACCGGCTAC
GGCTACCAGCAGGGCTACGGGCTGGCTATGGCGGCTACGACTACTCGCCCTATGGCTATTACGGCTACG
GCCCCGCTACGACTACAGTCAAGGTAGTACAACTACGGCAAGAGCCAGCGACGTGGTGCCATCAGAA
TAACTACAAGCCATAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC204360 representing NM_031266 Red=Cloning site Green=Tags(s) MSEAGEEQPMETTGTATENGHEAVPEGESPAAGTGAAAGAGGATAAPPSGNQNGAEGDQINASKNEEDAG KMFVGGLSWDTSKKDLKDYFTKFGVVDCTIKMDPNTGRSRGFGFILFKDAASVEKVLDQKEHRLDGRVI DPKKAMAMKKDPVKKIFVGGNLPEATEEKIREYFGEFGEIEAIELPMDPKLNKRRGFVFITFKEEEPVKK VLEKKFHTVSGSKCEIKVAQPKEVYQQQYGSGRGNRNRGNRSGGGGGGGGQSQSWNQGYGNYWNQGY GYQQGYGPGYGGYDYSYGYGYGPGYDYSQGSTNYGKSQRRGGHQNYPY TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_031266
ORF Size:	996 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_031266.2</u>
RefSeq Size:	1837 bp
RefSeq ORF:	999 bp
Locus ID:	3182
UniProt ID:	<u>Q99729</u>
Cytogenetics:	5q35.3
MW:	35.8 kDa