

Product datasheet for **RC204363L3V**

EVI5L (NM_145245) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	EVI5L (NM_145245) Human Tagged ORF Clone Lentiviral Particle
Symbol:	EVI5L
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_145245
ORF Size:	2382 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204363).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_145245.2 , NP_660288.1
RefSeq Size:	3834 bp
RefSeq ORF:	2385 bp
Locus ID:	115704
UniProt ID:	Q96CN4
Cytogenetics:	19p13.2
Domains:	TBC
Protein Families:	Druggable Genome
MW:	91.4 kDa


[View online »](#)

Gene Summary: Functions as a GTPase-activating protein (GAP) with a broad specificity.[UniProtKB/Swiss-Prot Function]