

Product datasheet for **RC204372A1V**

Human CALML5 (NM_017422) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human CALML5 (NM_017422) AAV Particle
Tag: Myc-DDK
Symbol: CALML5
Synonyms: CLSP
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC204372 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCCGGTGAGCTGACTCCTGAGGAGGAGGCCAGTACAAAAAGGCTTTCTCCGCGTTGACACGGATG
GAAACGGCACCATCAATGCCAGGAGCTGGGCGCGCGCTGAAGGCCACGGGCAAGAACCTCTCGGAGGC
CCAGCTAAGGAACTCATCTCCGAGGTTGACGGCGACGGCGACGGCGAAATCAGCTTCAGGAGTTCTCTG
ACGGCGGCAAGGAAGGCCAGGGCCGGCTGGAGGACCTGCAGGTCGCCTTCGCGCCTTCGACACGAGATG
GCGACGGCCACATCACCGTGGACGAGCTCAGGCGGGCCATGGCGGGGCTGGGGCAGCCGCTGCCGCAGGA
GGAGCTGGACGCCATGATCCGCGAGGCCGACGTGGACACGAGGACGGCGGGTGAACCTACGAGGAGTTCCGC
AGGATGCTCGCCAGGAG

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC204372 protein sequence
Red=Cloning site Green=Tags(s)

MAGELTPEEEAQYKAFSAVDTDGNGTINAQELGAALKATGKNLSEAQLRKLISEVDGDDGEISFQEF
TAARKARAGLEDLQVAFRAFDQDGDGHITVDELRRAMAGLGQPLPQEELDAMIREADVDDQGRVNYEEFA
RMLAQE

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2



[View online »](#)

ACCN:	NM_017422
ORF Size:	438 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_017422.3</u>
RefSeq Size:	893 bp
RefSeq ORF:	441 bp
Locus ID:	51806
UniProt ID:	<u>Q9NZT1</u>
Cytogenetics:	10p15.1
MW:	15.9 kDa