

Product datasheet for **RC204388L3V**

MTBP (NM_022045) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	MTBP (NM_022045) Human Tagged ORF Clone Lentiviral Particle
Symbol:	MTBP
Synonyms:	MDM2BP
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_022045
ORF Size:	2712 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204388).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_022045.2
RefSeq Size:	3087 bp
RefSeq ORF:	2715 bp
Locus ID:	27085
UniProt ID:	Q96DY7
Cytogenetics:	8q24.12
MW:	102.2 kDa



[View online »](#)

Gene Summary:

This gene encodes a protein that interacts with the oncoprotein mouse double minute 2. The encoded protein regulates progression through the cell cycle and may be involved in tumor formation. [provided by RefSeq, Aug 2012]