

Product datasheet for **RC204391L3V**

SLC2A6 (NM_017585) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	SLC2A6 (NM_017585) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SLC2A6
Synonyms:	GLUT6; GLUT9; HSA011372
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_017585
ORF Size:	1521 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204391).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_017585.2
RefSeq Size:	2563 bp
RefSeq ORF:	1524 bp
Locus ID:	11182
UniProt ID:	Q9UGQ3
Cytogenetics:	9q34.2
Protein Families:	Transmembrane
MW:	54.5 kDa



[View online »](#)

Gene Summary:

Hexose transport into mammalian cells is catalyzed by a family of membrane proteins, including SLC2A6, that contain 12 transmembrane domains and a number of critical conserved residues.[supplied by OMIM, Jul 2002]