

## Product datasheet for **RC204393L3V**

### **BXDC1 (RPF2) (NM\_032194) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	BXDC1 (RPF2) (NM_032194) Human Tagged ORF Clone Lentiviral Particle
Symbol:	BXDC1
Synonyms:	bA397G5.4; BXDC1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_032194
ORF Size:	918 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204393).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_032194.1</a>
RefSeq Size:	3704 bp
RefSeq ORF:	921 bp
Locus ID:	84154
UniProt ID:	<a href="#">Q9H7B2</a>
Cytogenetics:	6q21
MW:	35.6 kDa



[View online »](#)

**Gene Summary:**

Involved in ribosomal large subunit assembly. May regulate the localization of the 5S RNP/5S ribonucleoprotein particle to the nucleolus.[UniProtKB/Swiss-Prot Function]