

## Product datasheet for **RC204397L4V**

### **KIAA0859 (METTL13) (NM\_014955) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	KIAA0859 (METTL13) (NM_014955) Human Tagged ORF Clone Lentiviral Particle
Symbol:	KIAA0859
Synonyms:	5630401D24Rik; CGI-01; DFNB26; DFNB26M; DFNM1; feat; KIAA0859; METTL13
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_014955
ORF Size:	2100 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204397).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_014955.2</a> , <a href="#">NP_055770.1</a>
RefSeq Size:	3096 bp
RefSeq ORF:	1842 bp
Locus ID:	51603
UniProt ID:	<a href="#">Q8N6R0</a>
Cytogenetics:	1q24.3
Protein Families:	Druggable Genome
MW:	78.8 kDa



[View online »](#)

**Gene Summary:**

Dual methyltransferase that catalyzes methylation of elongation factor 1-alpha (EEF1A1 and EEF1A2) at two different positions, and is therefore involved in the regulation of mRNA translation (PubMed:30612740, PubMed:30143613). Via its C-terminus, methylates EEF1A1 and EEF1A2 at the N-terminal residue 'Gly-2' (PubMed:30143613). Via its N-terminus dimethylates EEF1A1 and EEF1A2 at residue 'Lys-55' (PubMed:30612740, PubMed:30143613). Has no activity towards core histones H2A, H2B, H3 and H4 (PubMed:30612740). [UniProtKB/Swiss-Prot Function]