

## Product datasheet for **RC204400L4V**

### POP5 (NM\_015918) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

|                                  |                                   |
|----------------------------------|-----------------------------------|
| <b>Product Type:</b>             | Lentiviral Particles              |
| <b>Symbol:</b>                   | POP5                              |
| <b>Synonyms:</b>                 | hPop5; HSPC004; RPP2; RPP20       |
| <b>Mammalian Cell Selection:</b> | Puromycin                         |
| <b>Vector:</b>                   | pLenti-C-mGFP-P2A-Puro (PS100093) |
| <b>Tag:</b>                      | mGFP                              |
| <b>ACCN:</b>                     | NM_015918                         |
| <b>ORF Size:</b>                 | 489 bp                            |

**ORF Nucleotide Sequence:** The ORF insert of this clone is exactly the same as(RC204400).

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

|                      |                             |
|----------------------|-----------------------------|
| <b>RefSeq:</b>       | <a href="#">NM_015918.3</a> |
| <b>RefSeq Size:</b>  | 801 bp                      |
| <b>RefSeq ORF:</b>   | 492 bp                      |
| <b>Locus ID:</b>     | 51367                       |
| <b>UniProt ID:</b>   | <a href="#">Q969H6</a>      |
| <b>Cytogenetics:</b> | 12q24.31                    |



**Protein Families:** Stem cell – Pluripotency

**MW:** 18.8 kDa

**Gene Summary:** Component of ribonuclease P, a protein complex that generates mature tRNA molecules by cleaving their 5'-ends (PubMed:11413139, PubMed:30454648). Also a component of the MRP ribonuclease complex, which cleaves pre-rRNA sequences (PubMed:28115465).  
[UniProtKB/Swiss-Prot Function]