

Product datasheet for **RC204411L3V**

FAIM1 (FAIM) (NM_001033032) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	FAIM1 (FAIM) (NM_001033032) Human Tagged ORF Clone Lentiviral Particle
Symbol:	FAIM1
Synonyms:	FAIM1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001033032
ORF Size:	537 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204411).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001033032.1
RefSeq Size:	1104 bp
RefSeq ORF:	540 bp
Locus ID:	55179
UniProt ID:	Q9NVQ4
Cytogenetics:	3q22.3
MW:	20.2 kDa



[View online »](#)

Gene Summary:

The protein encoded by this gene protects against death receptor-triggered apoptosis and regulates B-cell signaling and differentiation. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2011]