

Product datasheet for **RC204423A1V**

Human Kallikrein 7 (KLK7) (NM_005046) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human Kallikrein 7 (KLK7) (NM_005046) AAV Particle
Tag:	Myc-DDK
Symbol:	Kallikrein 7
Synonyms:	hK7; PRSS6; SCCE
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC204423 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCAAGATCCCTTCTCCTGCCCTGCAGATCCTACTGCTATCCTTAGCCTTGAAACTGCAGGAGAAG
AAGCCAGGGTGACAAGATTATTGATGGCGCCCCATGTGCAAGAGGCTCCACCCATGGCAGGTGGCCCT
GCTCAGTGGCAATCAGCTCCACTGCGGAGGCGTCTGGTCAATGAGCGCTGGGTGCTCACTGCCGCCAC
TGCAAGATGAATGAGTACACCGTGACCTGGGCAGTGATACGCTGGGCGACAGGAGAGCTCAGAGGATCA
AGGCCTCGAAGTCATTCGCCACCCCGGCTACTCCACACAGACCATGTTAATGACCTCATGCTCGTGAA
GCTCAATAGCCAGGCCAGGCTGTCATCCATGGTGAAGAAAGTCAGGCTGCCCTCCCGCTGCGAACCCCT
GGAACCACTGTACTGTCTCCGGCTGGGGCACTACCACGAGCCAGATGTGACCTTTCCCTCTGACCTCA
TGTGCGTGGATGTCAAGCTCATCTCCCCCAGGACTGCACGAAGGTTTACAAGGACTTACTGGAAAATTC
CATGCTGTGCGTGGCATCCCCGACTCCAAGAAAAACGCTGCAATGGTACTCAGGGGGACCGTTGGTG
TGCAGAGGTACCCTGCAAGGTCTGGTGTCTGGGGAACTTTCCCTTGGGGCAACCCAATGACCCAGGAG
TCTACACTCAAGTGTGAAGTTCACCAAGTGGATAAATGACACCATGAAAAGCATCGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC204423 protein sequence Red=Cloning site Green=Tags(s) MARSLLLPLQILLLSLALETAGEEAQGDKIIDGAPCARGSHPWQVALLSGNQLHCGGVLVNERWVLTAAH CKMNEYTVHLGSDTLGDRRAQRIKASKSFRHPGYSTQTHVNDLMLVKLNSQARLSSMVKKVRLPSRCEPP GTTCTVSGWGTTPSPDVTFFSDLMCVDKLISPQDCTKVYKDLLENSMLCAGIPDSKKNACNGDSGGPLV CRGTLQGLVSWGTFPWGQPNDPGVYTQVCKFTKWINDTMKKHR TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_005046
ORF Size:	759 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_005046.2</u>
RefSeq Size:	2104 bp
RefSeq ORF:	762 bp
Locus ID:	5650
UniProt ID:	<u>P49862</u>
Cytogenetics:	19q13.41
MW:	27.6 kDa