

Product datasheet for **RC204447L4V**

MICA (NM_000247) Human Tagged ORF Clone Lentiviral Particle

Product data:


Product Type:	Lentiviral Particles
Product Name:	MICA (NM_000247) Human Tagged ORF Clone Lentiviral Particle
Symbol:	MICA
Synonyms:	MIC-A; PERB11.1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_000247
ORF Size:	1149 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204447).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_000247.1
RefSeq Size:	1410 bp
RefSeq ORF:	1152 bp
Locus ID:	100507436
UniProt ID:	Q29983
Cytogenetics:	6p21.33
Domains:	MHC_I, IGc1
MW:	42.9 kDa



[View online »](#)

Gene Summary:

This gene encodes the highly polymorphic major histocompatibility complex class I chain-related protein A. The protein product is expressed on the cell surface, although unlike canonical class I molecules it does not seem to associate with beta-2-microglobulin. It is a ligand for the NKG2-D type II integral membrane protein receptor. The protein functions as a stress-induced antigen that is broadly recognized by intestinal epithelial gamma delta T cells. Variations in this gene have been associated with susceptibility to psoriasis 1 and psoriatic arthritis, and the shedding of MICA-related antibodies and ligands is involved in the progression from monoclonal gammopathy of undetermined significance to multiple myeloma. Alternative splicing of this gene results in multiple transcript variants. [provided by RefSeq, Jan 2014]

Product images:

[RC204447L4] was used to prepare Lentiviral particles using [TR30037] packaging kit. HEK293T cells were transduced with RC204447L4V particle to overexpress human MICA-mGFP fusion protein.