

Product datasheet for **RC204450L3V**

Cleavage stimulation factor 2 (CSTF2) (NM_001325) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Cleavage stimulation factor 2 (CSTF2) (NM_001325) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Cleavage stimulation factor 2
Synonyms:	CstF-64
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001325
ORF Size:	1731 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204450).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001325.2
RefSeq Size:	2027 bp
RefSeq ORF:	1734 bp
Locus ID:	1478
UniProt ID:	P33240
Cytogenetics:	Xq22.1
Domains:	RRM



[View online »](#)

MW: 61 kDa

Gene Summary: This gene encodes a nuclear protein with an RRM (RNA recognition motif) domain. The protein is a member of the cleavage stimulation factor (CSTF) complex that is involved in the 3' end cleavage and polyadenylation of pre-mRNAs. Specifically, this protein binds GU-rich elements within the 3'-untranslated region of mRNAs. [provided by RefSeq, Jul 2008]