

Product datasheet for **RC204456L4V**

DNA Polymerase gamma (POLG) (NM_002693) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	DNA Polymerase gamma (POLG) (NM_002693) Human Tagged ORF Clone Lentiviral Particle
Symbol:	DNA Polymerase gamma
Synonyms:	MDP1; MIRAS; MTDPS4A; MTDPS4B; PEO; POLG1; POLGA; SANDO; SCAE
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_002693
ORF Size:	3717 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204456).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_002693.1
RefSeq Size:	4464 bp
RefSeq ORF:	3720 bp
Locus ID:	5428
UniProt ID:	P54098
Cytogenetics:	15q26.1
Domains:	DNA_pol_A
Protein Families:	Druggable Genome



[View online »](#)

Protein Pathways: Metabolic pathways

MW: 139.6 kDa

Gene Summary: Mitochondrial DNA polymerase is heterotrimeric, consisting of a homodimer of accessory subunits plus a catalytic subunit. The protein encoded by this gene is the catalytic subunit of mitochondrial DNA polymerase. The encoded protein contains a polyglutamine tract near its N-terminus that may be polymorphic. Defects in this gene are a cause of progressive external ophthalmoplegia with mitochondrial DNA deletions 1 (PEOA1), sensory ataxic neuropathy dysarthria and ophthalmoparesis (SANDO), Alpers-Huttenlocher syndrome (AHS), and mitochondrial neurogastrointestinal encephalopathy syndrome (MNGIE). Two transcript variants encoding the same protein have been found for this gene. [provided by RefSeq, Jul 2008]