

### **Product datasheet for RC204461L1**

# FBP1 (NM\_000507) Human Tagged Lenti ORF Clone

#### **Product data:**

**Product Type:** Expression Plasmids

**Product Name:** FBP1 (NM\_000507) Human Tagged Lenti ORF Clone

Tag: Myc-DDK

Symbol: FBP1

Synonyms: FBP

Mammalian Cell None

Selection:

Vector:pLenti-C-Myc-DDK (PS100064)E. coli Selection:Chloramphenicol (34 ug/mL)

ORF Nucleotide The ORF insert of this clone is exactly the same as(RC204461).

Sequence:

**Restriction Sites:** Sgfl-Mlul

**Cloning Scheme:** 





<sup>\*</sup> The last codon before the Stop codon of the ORF.

**ACCN:** NM\_000507

ORF Size: 1014 bp



**OriGene Technologies, Inc.** 9620 Medical Center Drive, Ste 200

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com

CN: techsupport@origene.cn

### FBP1 (NM\_000507) Human Tagged Lenti ORF Clone - RC204461L1

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:** 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

**RefSeq:** <u>NM 000507.2</u>

RefSeq Size: 1647 bp
RefSeq ORF: 1017 bp
Locus ID: 2203
UniProt ID: P09467

Cytogenetics: 9q22.32

**Domains:** FBPase

**Protein Families:** Druggable Genome, Stem cell - Pluripotency

**Protein Pathways:** Fructose and mannose metabolism, Glycolysis / Gluconeogenesis, Insulin signaling pathway,

Metabolic pathways, Pentose phosphate pathway

MW: 36.8 kDa


**Gene Summary:** Fructose-1,6-bisphosphatase 1, a gluconeogenesis regulatory enzyme, catalyzes the

hydrolysis of fructose 1,6-bisphosphate to fructose 6-phosphate and inorganic phosphate. Fructose-1,6-diphosphatase deficiency is associated with hypoglycemia and metabolic


acidosis. [provided by RefSeq, Jul 2008]



# **Product images:**



Circular map for RC204461L1



Double digestion of RC204461L1 using Sgfl and Mlul  $\,$