

## Product datasheet for **RC204467L3V**

### **DHX30 (NM\_138615) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

<b>Product Type:</b>	Lentiviral Particles
<b>Product Name:</b>	DHX30 (NM_138615) Human Tagged ORF Clone Lentiviral Particle
<b>Symbol:</b>	DHX30
<b>Synonyms:</b>	DDX30; NEDMIAL; RETCOR
<b>Mammalian Cell Selection:</b>	Puromycin
<b>Vector:</b>	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
<b>Tag:</b>	Myc-DDK
<b>ACCN:</b>	NM_138615
<b>ORF Size:</b>	3582 bp
<b>ORF Nucleotide Sequence:</b>	The ORF insert of this clone is exactly the same as(RC204467).
<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>RefSeq:</b>	<a href="#">NM_138615.1</a>
<b>RefSeq Size:</b>	4104 bp
<b>RefSeq ORF:</b>	3585 bp
<b>Locus ID:</b>	22907
<b>UniProt ID:</b>	<a href="#">Q7L2E3</a>
<b>Cytogenetics:</b>	3p21.31
<b>Domains:</b>	DSRM, DEAD, helicase_C, HA2
<b>MW:</b>	133.9 kDa



[View online »](#)

**Gene Summary:**

DEAD box proteins, characterized by the conserved motif Asp-Glu-Ala-Asp (DEAD), are putative RNA helicases. They are implicated in a number of cellular processes involving alteration of RNA secondary structure such as translation initiation, nuclear and mitochondrial splicing, and ribosome and spliceosome assembly. Based on their distribution patterns, some members of this DEAD box protein family are believed to be involved in embryogenesis, spermatogenesis, and cellular growth and division. The family member encoded by this gene is a mitochondrial nucleoid protein that associates with mitochondrial DNA. It has also been identified as a component of a transcriptional repressor complex that functions in retinal development, and it is required to optimize the function of the zinc-finger antiviral protein. Alternatively spliced transcript variants have been found for this gene. [provided by RefSeq, Feb 2013]