

## Product datasheet for **RC204474L1V**

### Pannexin 1 (PANX1) (NM\_015368) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Pannexin 1 (PANX1) (NM_015368) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Pannexin 1
Synonyms:	MRS1; OOMD7; PX1; UNQ2529
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_015368
ORF Size:	1278 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204474).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_015368.3</a> , <a href="#">NP_056183.2</a>
RefSeq Size:	2782 bp
RefSeq ORF:	1281 bp
Locus ID:	24145
UniProt ID:	<a href="#">Q96RD7</a>
Cytogenetics:	11q21
Protein Families:	Transmembrane
MW:	47.9 kDa



[View online »](#)

**Gene Summary:**

The protein encoded by this gene belongs to the innexin family. Innexin family members are the structural components of gap junctions. This protein and pannexin 2 are abundantly expressed in central nerve system (CNS) and are coexpressed in various neuronal populations. Studies in *Xenopus* oocytes suggest that this protein alone and in combination with pannexin 2 may form cell type-specific gap junctions with distinct properties. [provided by RefSeq, Jul 2008]