

Product datasheet for **RC204510A1V**

Human RAB37 (NM_001006638) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human RAB37 (NM_001006638) AAV Particle
Tag: Myc-DDK
Symbol: RAB37
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC204510 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGACGGGCACGCCAGGCGCGGTTGCCACCCGGGATGGCGAGGCCCGGAGCGCTCCCCGCCCTGCAGTC
CGAGCTACGACCTCACGGGCAAGGTGATGCTTCTGGGAGACACAGGCGTCGGCAAAACATGTTTCCTGAT
CCAATTCAAAGACGGGCGCTTCTGTCCGGAACCTTCATAGCCACCGTCGGCATAGACTTCAGGAACAAG
GTGGTGACTGTGGATGGCGTGAGAGTGAAGCTGCAGATCTGGGACACCGCTGGGAGGAACGGTTCGAA
GCGTCACCCATGCTTATTACAGAGATGCTCAGGCCTTGCTTCTGCTGTATGACATCACCACAAATCTTC
TTTCGACAACATCAGGGCCTGGCTCACTGAGATTCATGAGTATGCCAGAGGGACGTGGTGATCATGCTG
CTAGGCAACAAGGCGGATATGAGCAGCGAAAGAGTGATCCGTTCCGAAGACGGAGAGACCTTGGCCAGGG
AGTACGGTGTTCCCTTCTGGAGACCAGCGCCAAGACTGGCATGAATGTGGAGTTAGCCTTTCTGGCCAT
CGCCAAGGAAGTAAATACCGGGCCGGGCATCAGGCGGATGAGCCAGCTTCAGATCCGAGACTATGTA
GAGTCCCAGAAGAAGCGCTCCAGCTGCTGCTCCTTCATG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC204510 protein sequence
Red=Cloning site Green=Tags(s)

MTGTPGAVATRDGEAPERSPPCSPSYDLTGKVMLLGDTGVGKTCFLIQFKDGAFLSGTFIATVGIDFRNK
VVTVDGVRVKLIWDTAGQERFRSVTHAYYRDAQALLLYDITNKSSFNDNIRAWLTEIHEYAQRDVVIML
LGNAKDMSSSERVIRSEDGETLAREYGVFPLETSAKTGMNVELAFLAIKELKYRAGHQADEPSFQIRDYV
ESQKKRSSCCSFM

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human



[View online »](#)

Serotype:	AAV-2
ACCN:	NM_001006638
ORF Size:	669 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001006638.1</u>
RefSeq Size:	2698 bp
RefSeq ORF:	672 bp
Locus ID:	326624
UniProt ID:	<u>Q96AX2</u>
Cytogenetics:	17q25.1
MW:	24.8 kDa