

Product datasheet for **RC204518L4V**

Epac1 (RAPGEF3) (NM_006105) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Epac1 (RAPGEF3) (NM_006105) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Epac1
Synonyms:	bcm910; CAMP-GEFI; EPAC; EPAC1; HSU79275
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_006105
ORF Size:	2772 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204518).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_006105.3 , NP_006096.2
RefSeq Size:	5773 bp
RefSeq ORF:	2646 bp
Locus ID:	10411
UniProt ID:	O95398
Cytogenetics:	12q13.11
Domains:	DEP, cNMP, RasGEFN, RasGEF
Protein Pathways:	Leukocyte transendothelial migration, Long-term potentiation



[View online »](#)

MW: 103.7 kDa

Gene Summary: Guanine nucleotide exchange factor (GEF) for RAP1A and RAP2A small GTPases that is activated by binding cAMP. Through simultaneous binding of PDE3B to RAPGEF3 and PIK3R6 is assembled in a signaling complex in which it activates the PI3K gamma complex and which is involved in angiogenesis. Plays a role in the modulation of the cAMP-induced dynamic control of endothelial barrier function through a pathway that is independent on Rho-mediated signaling. Required for the actin rearrangement at cell-cell junctions, such as stress fibers and junctional actin.[UniProtKB/Swiss-Prot Function]