

## Product datasheet for **RC204528L3V**

### **RNF113A (NM\_006978) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	RNF113A (NM_006978) Human Tagged ORF Clone Lentiviral Particle
Symbol:	RNF113A
Synonyms:	Cwc24; RNF113; TTD5; ZNF183
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_006978
ORF Size:	1029 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204528).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_006978.1</a>
RefSeq Size:	1349 bp
RefSeq ORF:	1032 bp
Locus ID:	7737
UniProt ID:	<a href="#">O15541</a>
Cytogenetics:	Xq24
Domains:	zf-CCCH, RING
Protein Families:	Druggable Genome, Protease



[View online »](#)

**MW:** 38.6 kDa

**Gene Summary:** This intronless gene encodes a protein which contains a C3H1-type zinc finger domain and a C3HC4 Ring-type (Really Interesting New Gene-type) zinc finger domain. The Ring-type zinc finger domain is identified in various tumor suppressors, DNA repair genes and cytokine receptor-associated molecules, and is probably involved in mediating protein-protein interactions. [provided by RefSeq, May 2010]