

## Product datasheet for **RC204545L2V**

### CYC1 (NM\_001916) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	CYC1 (NM_001916) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CYC1
Synonyms:	MC3DN6; UQCR4
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_001916
ORF Size:	975 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204545).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001916.2</a> , <a href="#">NP_001907.2</a>
RefSeq Size:	1251 bp
RefSeq ORF:	978 bp
Locus ID:	1537
UniProt ID:	<a href="#">P08574</a>
Cytogenetics:	8q24.3
Domains:	Cytochrome_C1
Protein Families:	Druggable Genome



[View online »](#)

**Protein Pathways:** Alzheimer's disease, Cardiac muscle contraction, Huntington's disease, Metabolic pathways, Oxidative phosphorylation, Parkinson's disease

**MW:** 35.4 kDa

**Gene Summary:** This gene encodes a subunit of the cytochrome bc1 complex, which plays an important role in the mitochondrial respiratory chain by transferring electrons from the Rieske iron-sulfur protein to cytochrome c. Mutations in this gene may cause mitochondrial complex III deficiency, nuclear type 6. [provided by RefSeq, Dec 2013]