

Product datasheet for RC204548L1

GFAP (NM_002055) Human Tagged Lenti ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: GFAP (NM_002055) Human Tagged Lenti ORF Clone

Tag: Myc-DDK

Symbol: GFAP

Synonyms: ALXDRD

Mammalian Cell None

Selection:

Vector:pLenti-C-Myc-DDK (PS100064)E. coli Selection:Chloramphenicol (34 ug/mL)

ORF Nucleotide The ORF insert of this clone is exactly the same as(RC204548).

Sequence:

Restriction Sites: Sgfl-Mlul

Cloning Scheme:





st The last codon before the Stop codon of the ORF.

ACCN: NM_002055

ORF Size: 1296 bp



OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com



GFAP (NM_002055) Human Tagged Lenti ORF Clone - RC204548L1

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 002055.2</u>

 RefSeq Size:
 3097 bp

 RefSeq ORF:
 1299 bp

 Locus ID:
 2670

 UniProt ID:
 P14136

Cytogenetics: 17q21.31

Domains: filament, filament_head

Protein Families: ES Cell Differentiation/IPS

MW: 49.9 kDa


Gene Summary: This gene encodes one of the major intermediate filament proteins of mature astrocytes. It is

used as a marker to distinguish astrocytes from other glial cells during development. Mutations in this gene cause Alexander disease, a rare disorder of astrocytes in the central nervous system. Alternative splicing results in multiple transcript variants encoding distinct


isoforms. [provided by RefSeq, Oct 2008]



Product images:



Circular map for RC204548L1



Double digestion of RC204548L1 using Sgfl and Mlul $\,$